



AN323 用加速溶剂萃取 (ASE) 技术提取环境样品中的多氯二苯二噁英和多氯二苯呋喃

PCDDs: 多氯二苯二噁英

OCPs: 有机氯农药

PCDFs: 多氯二苯呋喃

OPPS: 有机磷农药

样品基体: 环境样品 (烟囱砖灰, 城市尘埃, 飞尘和沉积物)

仪器: Dionex ASE200

GC/MS

低压 LC (样品清洗)

溶剂: 甲苯 (农药级), 乙酸 (试剂级), HCl (试剂级)

推荐的萃取条件

压力: 13.8Mpa (2000psi)

温度: 175—200°C

静态时间: 5—15min

溶剂: 甲苯, 或甲苯/乙酸 (5%, V/V) (如果用 HCl 作前处理)

冲洗体积: 60—70%

N₂ 吹扫时间: 60—100 秒

静态循环次数: 2 或 3 次

分析结果:

Table 2 Average Values (ng/kg) from Ground Chimney Brick—Comparison of Soxhlet vs. ASE		
Group Totals	Soxhlet (n=1)	ASE* (n=2)
Total T ₁ CDD	440	530
Total P ₁ CDD	900	940
Total H ₁ CDD	1800	2000
Total H ₂ CDD	2000	2100
Total O ₁ CDD	2900	2600
Total T ₁ CDF	2300	2600
Total P ₁ CDF	4100	4300
Total H ₁ CDF	4700	4700
Total H ₂ CDF	2800	2600
Total O ₁ CDF	2000	2000
Congeners	Soxhlet (n=1)	ASE* (n=2)
2,3,7,8-T ₁ CDD	6.0	6.0
1,2,3,7,8-P ₁ CDD	52	57
1,2,3,4,7,8-H ₁ CDD	46	52
1,2,3,6,7,8-H ₂ CDD	120	130
1,2,3,7,8,9-H ₂ CDD	97	1000
1,2,3,4,6,7,8-H ₃ CDD	1000	820
2,3,7,8-T ₁ CDF	160	180
1,2,3,7,8-(+1,2,3,4,8)-P ₁ CDF	430	470
2,3,4,7,8-P ₁ CDF	390	390
1,2,3,4,7,8-(+1,2,3,4,7,9)-H ₂ CDF	1100	1100
1,2,3,6,7,8-H ₂ CDF	540	570
2,3,4,6,7,8-H ₃ CDF	400	360
1,2,3,7,8,9-H ₂ CDF	42	42
1,2,3,4,6,7,8-H ₃ CDF	2100	2000
1,2,3,4,7,8,9-H ₃ CDF	140	120
Toxicity equivalent (NATO)	540	540
Toxicity equivalent (BgW)	490	510

Values are corrected for recovery of ¹⁴C-labeled surrogates.
*Sum of two extractions of each sample.

Table 3 Average Values (ng/kg) from Urban Dust—Comparison of Soxhlet vs. ASE		
Group Totals	Soxhlet (n=1)	ASE* (n=2)
Total T ₁ CDD	182	325
Total P ₁ CDD	175	221
Total H ₁ CDD	86.7	81.7
Total H ₂ CDD	221	217
Total O ₁ CDD	445	314
Total T ₁ CDF	333	419
Total P ₁ CDF	146	179
Total H ₁ CDF	65.9	122
Total H ₂ CDF	13.2	29.4
Total O ₁ CDF	n.d. (10)	n.d. (10)
Congeners	Soxhlet (n=1)	ASE* (n=2)
2,3,7,8-T ₁ CDD	3.3	3.2
1,2,3,7,8-P ₁ CDD	11.8	13.1
1,2,3,4,7,8-H ₁ CDD	9.8	8.0
1,2,3,6,7,8-H ₂ CDD	11.5	9.5
1,2,3,7,8,9-H ₂ CDD	n.d. (8)	n.d. (8)
1,2,3,4,6,7,8-H ₃ CDD	113	107
2,3,7,8-T ₁ CDF	12.5	18.6
1,2,3,7,8-(+1,2,3,4,8)-P ₁ CDF	9.9	12.0
2,3,4,7,8-P ₁ CDF	13.9	18.1
1,2,3,4,7,8-(+1,2,3,4,7,9)-H ₂ CDF	18.7	23.7
1,2,3,6,7,8-H ₂ CDF	10.7	15.8
2,3,4,6,7,8-H ₃ CDF	3.3	8.7
1,2,3,7,8,9-H ₂ CDF	n.d. (2)	n.d. (2)
1,2,3,4,6,7,8-H ₃ CDF	13.2	29.4
1,2,3,4,7,8,9-H ₃ CDF	n.d. (3)	n.d. (3)
Toxicity equivalent (NATO)	25.0	28.7
Toxicity equivalent (BgW)	24.2	29.2

Values are corrected for recovery of ¹⁴C-labeled surrogates.
n.d.—not detected. Detection limit, in ppt, given in parentheses.
*Sum of two extractions of each sample.

**Table 4 Average Values (µg/kg) from Fly Ash^a
Comparison of Soxhlet vs. ASE (HCl Pretreatment)**

Group Totals	Soxhlet (n=1)	ASE ^b (n=2)
Total T ₄ CDD	12.0	12.4
Total P ₅ CDD	16.6	20.5
Total H ₆ CDD	38.2	42.4
Total H ₇ CDD	15.0	19.8
Total O ₈ CDD	11.4	12.8
Total T ₄ CDF	60.5	67.5
Total P ₅ CDF	83.5	87.3
Total H ₆ CDF	65.2	73.5
Total H ₇ CDF	28.1	32.2
Total O ₈ CDF	13.5	15.8
Congeners	Soxhlet (n=1)	ASE ^b (n=2)
2,3,7,8-T ₄ CDD	0.32	0.36
1,2,3,7,8-P ₅ CDD	1.6	2.1
1,2,3,4,7,8-H ₆ CDD	1.2	1.4
1,2,3,6,7,8-H ₆ CDD	2.4	2.7
1,2,3,7,8,9-H ₆ CDD	2.4	2.3
1,2,3,4,6,7,8-H ₇ CDD	8.2	9.6
2,3,7,8-T ₄ CDF	3.7	4.3
1,2,3,7,8(+1,2,3,4,8)-P ₅ CDF	4.2	4.6
2,3,4,7,8-P ₅ CDF	5.6	6.6
1,2,3,4,7,8(+1,2,3,4,7,9)-H ₆ CDF	7.8	8.7
1,2,3,6,7,8-H ₆ CDF	7.2	8.5
2,3,4,6,7,8-H ₆ CDF	6.6	7.2
1,2,3,7,8,9-H ₆ CDF	0.43	0.56
1,2,3,4,6,7,8-H ₇ CDF	18.0	17.6
1,2,3,4,7,8,9-H ₇ CDF	2.3	2.4
Toxicity equivalent (NATO)	7.6	8.9
Toxicity equivalent (BgVV)	7.3	14.2

Values are corrected for recovery of ¹⁴C-labeled surrogates.^aPretreatment with HCl for 30 min followed by water rinse.^bSum of two extractions of each sample.**Table 5 Average Values (µg/kg) from Fly Ash^a
Comparison of Soxhlet (HCl Pretreatment) vs. ASE
(No HCl Pretreatment; Toluene/Acetic Acid Solvent)**

Group Totals	Soxhlet (n=1)	ASE ^b (n=2)
Total T ₄ CDD	12.0	10.5
Total P ₅ CDD	16.6	16.2
Total H ₆ CDD	38.2	36.7
Total H ₇ CDD	15.0	16.0
Total O ₈ CDD	11.4	10.6
Total T ₄ CDF	60.5	56.1
Total P ₅ CDF	83.5	77.4
Total H ₆ CDF	65.2	46.1
Total H ₇ CDF	28.1	26.5
Total O ₈ CDF	13.5	13.9
Congeners	Soxhlet (n=1)	ASE ^b (n=2)
2,3,7,8-T ₄ CDD	0.32	0.28
1,2,3,7,8-P ₅ CDD	1.6	1.7
1,2,3,4,7,8-H ₆ CDD	1.2	1.2
1,2,3,6,7,8-H ₆ CDD	2.4	2.4
1,2,3,7,8,9-H ₆ CDD	2.4	2.2
1,2,3,4,6,7,8-H ₇ CDD	8.2	8.1
2,3,7,8-T ₄ CDF	3.7	3.4
1,2,3,7,8(+1,2,3,4,8)-P ₅ CDF	4.2	3.9
2,3,4,7,8-P ₅ CDF	5.6	5.8
1,2,3,4,7,8(+1,2,3,4,7,9)-H ₆ CDF	7.8	5.4
1,2,3,6,7,8-H ₆ CDF	7.2	5.3
2,3,4,6,7,8-H ₆ CDF	6.6	4.5
1,2,3,7,8,9-H ₆ CDF	0.43	0.30
1,2,3,4,6,7,8-H ₇ CDF	18.0	16.8
1,2,3,4,7,8,9-H ₇ CDF	2.3	2.0
Toxicity equivalent (NATO)	7.6	7.0
Toxicity equivalent (BgVV)	7.3	6.3

Values are corrected for recovery of ¹⁴C-labeled surrogates.^aPretreatment with HCl for 30 min followed by water rinse for Soxhlet samples.^bNo HCl pretreatment, and toluene/acetic acid solvent used for ASE samples.^cSum of two extractions of each sample.

**Table 6 Average Values (ng/kg) from EC-2—Comparison of Soxhlet vs. ASE**

Group Totals	Soxhlet Extraction (n=10)			Accelerated Solvent Extraction (n=2)		
	Value	% RSD	Isomers	Value	% RSD	Isomers
Total T ₄ CDD	430	9.7	8	370	1.9	9
Total P ₅ CDD	300	3.7	11	280	7.7	11
Total H ₆ CDD	720	5.8	7	690	2.0	7
Total H ₇ CDD	1300	7.0	2	1300	0.0	2
Total O ₈ CDD	4000	6.2	1	4200	0.0	1
Total T ₄ CDF	620	12	17	380	19	19
Total P ₅ CDF	820	9.4	14	710	7.0	17
Total H ₆ CDF	1900	5.7	12	1900	0.0	13
Total H ₇ CDF	3800	8.2	4	3900	3.6	4
Total O ₈ CDF	7800	8.3	1	7000	3.1	1
Congeners	Value	% RSD	% Recovery ^a	Value	% RSD	% Recovery ^a
2,3,7,8-T ₄ CDD	270	9.1	68	270	0.0	72
1,2,3,7,8-P ₅ CDD	24	12	74	22	3.3	81
1,2,3,4,7,8-H ₆ CDD	23	8.3	76	24	3.0	80
1,2,3,6,7,8-H ₆ CDD	83	3.6	78	87	0.8	54
1,2,3,7,8,9-H ₆ CDD	60	6.2	77	57	7.4	67
1,2,3,4,6,7,8-H ₇ CDD	720	6.7	81	720	1.0	79
2,3,7,8-T ₄ CDF ^b	100	7.3	68	82	2.6	70
1,2,3,7,8-P ₅ CDF	39	14	74	36	3.9	76
2,3,4,7,8-P ₅ CDF	62	5.5	79	60	0.0	75
1,2,3,4,7,8-H ₆ CDF	740	5.3	81	690	0.0	70
1,2,3,6,7,8-H ₆ CDF	120	6.2	81	120	0.0	50
2,3,4,6,7,8-H ₆ CDF	45	9.0	82	60	1.2	69
1,2,3,7,8,9-H ₆ CDF	4.9	31	84	5.3	15	70
1,2,3,4,6,7,8-H ₇ CDF	2600	6.7	85	2500	0.0	74
1,2,3,4,7,8,9-H ₇ CDF	160	5.5	83	160	0.0	72

Values are corrected for recovery of ¹⁴C-labeled surrogates.^aRefers to recovery of corresponding ¹⁴C-labeled surrogate.^bMaximum possible concentration due to potential chromatographic overlap.

**Table 7 Average Values (ng/kg) from HS-2—Comparison of Soxhlet vs. ASE**

Group Totals	Soxhlet Extraction (n=4)			Accelerated Solvent Extraction (n=2)		
	Value	% RSD	Isomers	Value	% RSD	Isomers
Total T ₄ CDD	3.9	14	2	2.5	34	5
Total P ₅ CDD	17	7.8	6	10	10	9
Total H ₆ CDD	510	5.6	8	570	1.3	7
Total H ₇ CDD	4700	8.3	2	5100	11	2
Total O ₈ CDD	6500	4.2	1	7100	0.0	1
Total T ₄ CDF	39	11	13	24	3.0	14
Total P ₅ CDF	33	13	8	28	0.0	11
Total H ₆ CDF	89	3.2	6	87	12	10
Total H ₇ CDF	293	3.3	4	310	0.0	4
Total O ₈ CDF	300	3.8	1	280	2.6	1
Congeners	Value	% RSD	% Recovery ^a	Value	% RSD	% Recovery ^a
2,3,7,8-T ₄ CDD	n.d. (1)		62	n.d. (1)		71
1,2,3,7,8-P ₅ CDD	1.6	4.6	69	n.d. (1)		75
1,2,3,4,7,8-H ₆ CDD	4.5	4.8	74	5.2	11	73
1,2,3,6,7,8-H ₆ CDD	19	4.3	75	21	0.0	50
1,2,3,7,8,9-H ₆ CDD	24	4.3	74	28	2.6	61
1,2,3,4,6,7,8-H ₇ CDD	1200	8.1	80	1300	0.0	93
2,3,7,8-T ₄ CDF ^b	8.5	11	62	6.6	5.4	65
1,2,3,7,8-P ₅ CDF	1.9	17	68	2.0	0.0	72
2,3,4,7,8-P ₅ CDF	3.7	7.9	71	3.7	3.8	59
1,2,3,4,7,8-H ₆ CDF	17	7.3	79	17	4.3	70
1,2,3,6,7,8-H ₆ CDF	3.7	5.6	80	4.0	5.4	49
2,3,4,6,7,8-H ₆ CDF	3.7	18	81	4.4	3.2	61
1,2,3,7,8,9-H ₆ CDF	n.d. (1)		83	n.d. (1)	0.0	75
1,2,3,4,6,7,8-H ₇ CDF	91	1.6	83	96	3.7	82
1,2,3,4,7,8,9-H ₇ CDF	5.2	6.7	84	5.3	6.7	83

n.d.—not detected. Detection limit, in ppt, given in parentheses. Values not used for statistical calculations.

Values are corrected for recovery of ¹⁴C-labeled surrogates.^aRefers to recovery of corresponding ¹⁴C-labeled surrogates.^bMaximum possible concentration due to potential chromatographic overlap.

**Table 8 Average Values (ng/kg) from Highly Contaminated Sources—Comparison of Soxhlet vs. ASE**

Group Totals	Hamilton Harbor		Parrots Bay	
	Soxhlet	ASE	Soxhlet	ASE
Total T ₄ CDD	50 ⁹	14 ⁹	39 ⁹	48 ⁹
Total P ₅ CDD	63 ¹²	15 ⁹	87 ¹⁰	66 ¹⁰
Total H ₆ CDD	220 ⁷	180 ⁷	230 ⁶	200 ⁷
Total H ₇ CDD	850 ²	810 ²	580 ²	530 ²
Total O ₈ CDD	3100	3100	1900	1600
Total T ₄ CDF	370 ¹⁷	130 ^{11a}	400 ¹⁴	270 ¹¹
Total P ₅ CDF	290 ¹³	110 ^{11a}	180 ⁸	170 ¹¹
Total H ₆ CDF	240 ¹⁴	160 ¹¹	230 ⁶	230 ¹⁰
Total H ₇ CDF	350 ⁴	290 ⁴	400 ⁴	360 ⁴
Total O ₈ CDF	270	210	510	370
Congeners				
2,3,7,8-T ₄ CDD	3.7	3.1	19	19
1,2,3,7,8-P ₅ CDD	5.1	5.4	8.3	6.0
1,2,3,4,7,8-H ₆ CDD	6.4	7.2	8.6	6.7
1,2,3,6,7,8-H ₆ CDD	27	26	26	17
1,2,3,7,8,9-H ₆ CDD	20	28	24	18
1,2,3,4,6,7,8-H ₇ CDD	460	430	280	250
2,3,7,8-T ₄ CDF ^a	61	44 ^a	80	48
1,2,3,7,8-P ₅ CDF	14	14	n.d. (20)	9.8
2,3,4,7,8-P ₅ CDF	26	25 ^a	22	14
1,2,3,4,7,8-H ₆ CDF	27	37	79	59
1,2,3,6,7,8-H ₆ CDF	17	16	n.d. (20)	15
2,3,4,6,7,8-H ₆ CDF	14	14	21	11
1,2,3,7,8,9-H ₆ CDF	n.d. (2)	1.6	4.9	n.d. (1)
1,2,3,4,6,7,8-H ₇ CDF	130	130	270	220
1,2,3,4,7,8,9-H ₇ CDF	14	13	17	12

Values are corrected for recovery of isotopically labeled surrogate standards.

n.d.—not detected. Detection limit, in ppt, given in parentheses.

Superscripts indicate the number of isomers detected.

^aRecoveries outside the range 25% to 150%; results are not to be used for regulatory compliance purposes.¹⁴Maximum possible concentration due to potential chromatographic overlap.